# **Hydrogen Generators**

For GC & GC/MS fuel and carrier gas applications

## Hydrogen on demand up to 1,300 ml/min

The Parker Balston ultra high purity hydrogen generators offer the optimum combination of safe operation, reliability, performance and low cost of ownership.

Utilising field proven PEM cell technology, hydrogen is produced on demand from deionized water and electricity. It then passes through a palladium purifier module to remove oxygen down to <0.01ppm. Innovative control software allows remote monitoring and the connection of several generators together (cascading) to supply applications with a large hydrogen demand.

With an output capacity of up to 1,300 ml/min at a purity of 99.99999+% and pressures up to 11.9 bar, the generator is ideal for all GC & GC/MS carrier gas (including fast GC) and combustion detector applications, and a more cost effective alternative to helium. To ensure continuous operation all generators include automatic water feed as standard.

## **Contact Information:**

Parker Hannifin Manufacturing Limited Filtration and Separation Division Dukesway, Team Valley Trading Estate Gateshead, Tyne and Wear England NE11 0PZ

Tel: +44 (0)191 402 9000 Fax: +44 (0)191 482 6296 Email: balstonukinfo@parker.com www.parker.com/dhFNS



## **Product Features:**

- Flow capacities up to 1,300 ml/min and pressures up to 11.9 bar
- Ultra high purity, 99.99999+%, an ideal alternative to helium for GC & GC/MS carrier gas applications
- 3 year standard cell warranty
- Designed to run continuously with automatic water feed as standard
- Automatic shut-down in case of a hydrogen leak
- Minimal maintenance and operator attention





### **Certified safety**

Parker Balston hydrogen generators utilise an exclusive Proton Exchange Membrane and palladium catalyst to produce ultra high purity hydrogen on demand.

A built in sensing circuit shuts the generator down if a hydrogen leak is detected and an eight stage explosion protection systems ensures the highest level of operator safety.

A sophisticated control system connected to a liquid crystal display continuously monitors the vital operating parameters to ensure a safe and consistent performance.

That's why the Parker Balston hydrogen generators meet the strict safety guidelines to be certified for CE, CSA and UL approval.

#### **Proven technology**

Parker Balston's exclusive Proton Exchange Membrane is proven in thousands of GC installations worldwide. Maintenance requires only a few moments per year, so no inconvenient extended downtime, simply change the deionizer cartridge and filter every 6 months.

High pressure for fast GC is achieved via an electrolytic cell utilising patented interlocking gaskets. The cell is protected by alarms that will detect low water quality and low water levels and a small internal pump ensures the cell is supplied with water at all times.

Parker Balston's hydrogen generators are the most reliable hydrogen generators on the market today.

| Prin | cipal | Spec | ifica | tion |
|------|-------|------|-------|------|
|------|-------|------|-------|------|

| Model                          | H2PEMPD-510          | H2PEMPD-650          | H2PEMPD-850          | H2PEMPD-1100         | H2PEMPD-1300             |
|--------------------------------|----------------------|----------------------|----------------------|----------------------|--------------------------|
| Purity                         | 99.99999+%           | 99.99999+%           | 99.99999+%           | 99.99999+%           | 99.99999+%               |
| Flow Rates                     | 510 ml/min           | 650 ml/min           | 850 ml/min           | 1,100 ml/min         | 1,300 ml/min             |
| Outlet Connection              | 1/4" compression         |
| Delivery Pressure (adjustable) | 11.9                 | 11.9                 | 11.9                 | 11.9                 | 6.8 or 11.9 bar options* |
| Remote Monitoring              | Yes                  | Yes                  | Yes                  | Yes                  | Yes                      |
| Auto Water Fill                | Included as standard     |
| Water Quality Required         | >5 Mohm                  |
| Ambient Temperature            | 10 to 35°C               |
| Electrical Requirements        | 100-230V - 50/60Hz       |
| Power Consumption              | 270 Watts            | 300 Watts            | 340 Watts            | 420 Watts            | 450 Watts                |
| Dimensions (H x W x D)         | 435 x 342 x 533 mm   |                      |                      |                      |                          |
| Weight (Shipping)              | 25 kg (29)               |

\* for 6.8 bar option add suffix -100 and 11.9 bar option suffix -175 ie H2PEMPD-850-100 - see ordering details below.

#### **Ordering Information**

| Description                                | Model Number     | Maintenance Items   | Model Number  | Change Frequency |  |
|--|------------------|---|---------------|------------------|--|
| 1,300 ml/min Hydrogen Generator - 6.8 bar  | H2PEMPD-1300-100 | 6 Month Maintenance Kit   | MKH2PEMPD-6M  | 6 Months         |  |
| 510 ml/min Hydrogen Generator - 11.9 bar   | H2PEMPD-510-175  | 24 Month Maintenance Kit  | MKH2PEMPD-24M | 24 Months        |  |
| 650 ml/min Hydrogen Generator - 11.9 bar   | H2PEMPD-650-100  |   |               |                  |  |
| 850 ml/min Hydrogen Generator - 11.9 ba    | H2PEMPD-850-175  |   |               |                  |  |
| 1,100 ml/min Hydrogen Generator - 11.9 bar | H2PEMPD-1100-175 |   |               |                  |  |
| 1,300 ml/min Hydrogen Generator - 11.9 bar | H2PEMPD-1300-175 |   |               |                  |  |
|  |                  |   |               |                  |  |
|  |                  |   |               |                  |  |
|  |                  |   |               |                  |  |
|  |                  |   |               |                  |  |
|  | B02-0795         |   |               |                  |  |
| H2PEMPD - Remote Interface                 |                  |   |               |                  |  |
| H2PEMPD - 8 Port Manifold                  | B02-0853         |   |               |                  |  |
| H2PEMPD - H2 Load Balance Assy             | B02-0851         | All LOREMOD serverses are filled with oute water fill (AMC) as standard Miles water   |               |                  |  |
| H2PEMPD - Kit, Hydrogen Purge              | B02-0852         | All H2PEMPD generators are fitted with auto water fill (AWF) as standard. When using<br>the AWF the maximum water inlet pressure is 0.1 bar. The B02-0880 auto water feed |               |                  |  |
| H2PEMPD - 2 Port Manifold                  | B02-0861         | regulator assemble is available if required. The B02-0880 is supplied as standard with the H2PEMPD-1100 and H2PEMPD-1300  |               |                  |  |

©2015 Parker Hannifin Corporation. All rights reserved.

Catalogue: S3.2.466\_EN 04/15



ENGINEERING YOUR SUCCESS.